

## Making 5' Probe for Conditional Knock-Out Mice

1. Digest Eco B/E subclone with Pst I and Xba I overnight at 37°C.

| Total volume: 200 $\mu$ L          |             | DNA conc.: 10 $\mu$ g / 50 $\mu$ l |                | Enzyme: 200 Units             |                       |
|------------------------------------|-------------|------------------------------------|----------------|-------------------------------|-----------------------|
| Subclone Conc. ( $\mu$ g/ $\mu$ l) | Enzyme      | Enzyme ( $\mu$ l)                  | DNA ( $\mu$ L) | ddH <sub>2</sub> O ( $\mu$ l) | 10x Buffer ( $\mu$ l) |
|                                    | Pst I, H.C. | 5                                  |                |                               | 20                    |
| --                                 | Xba I, H.C. | 5                                  | --             | --                            | --                    |

2. Run analytical gel to check digest.
  - a) Use 1% agarose gel with 5  $\mu$ l Ethidium Bromide / 100  $\mu$ l agarose gel.
  - b) Use mini-gel apparatus (~50 ml total gel volume).
  - c) Prepare samples by mixing together the following:
    - 2  $\mu$ l DNA sample
    - 1  $\mu$ l Orange G
    - 7  $\mu$ l ddH<sub>2</sub>O
  - d) Use 2  $\mu$ l 1 Kb Plus DNA ladder.
3. If digested completely, run samples on medium-sized gel (~100 mL total gel volume).
  - a) Use 1% agarose gel with 5  $\mu$ l Ethidium Bromide / 100  $\mu$ l agarose gel.
  - b) Add 1  $\mu$ l Orange G / 10  $\mu$ l digested sample.
  - c) Use 5  $\mu$ l 1 Kb Plus DNA ladder.
  - d) Run gel for ~5 hr at 100 V.
4. Excise 0.6 Kb band from gel (\_\_\_ band from bottom). Take pictures of before & after.
5. Extract DNA from Gel using Qiagen QIAEX II Agarose Gel Extraction Protocol.
6. Quantify sample using lambda ladder.
  - a) Use 1% agarose gel with 5  $\mu$ l Ethidium Bromide / 100  $\mu$ l agarose gel.
  - b) Use mini-gel apparatus (~50 ml total gel volume).
  - c) Prepare samples by mixing together the following:
    - 5  $\mu$ l DNA sample
    - 2  $\mu$ l Orange G
    - 3  $\mu$ l ddH<sub>2</sub>O
  - d) Prepare  $\lambda$  ladder by mixing together the following and incubating for 10 minutes at 65°C:
    - 1  $\mu$ l  $\lambda$  marker
    - 2  $\mu$ l Orange G
    - 7  $\mu$ l ddH<sub>2</sub>O
  - e) Use 5  $\mu$ l  $\lambda$  ladder.
7. Compare sample band with  $\lambda$  ladder bands to determine concentration.