

## Making 3' Probe for Knock-In Mice

1. Digest Bam H/B subclone with Xho I and Hind III overnight at 37°C.

Total volume: 200 $\mu$ L		DNA conc.: 10 $\mu$ g / 50 $\mu$ l		Enzyme: 200 Units	
Subclone Conc. ( $\mu$ g/ $\mu$ l)	Enzyme	Enzyme ( $\mu$ l)	DNA ( $\mu$ L)	ddH <sub>2</sub> O ( $\mu$ l)	10x Buffer ( $\mu$ l)
	Hind III, H.C.	5			20
--	Xho I, H.C.	5	--	--	--

2. Run analytical gel to check digest.
  - a) Use 1% agarose gel with 5  $\mu$ l Ethidium Bromide / 100  $\mu$ l agarose gel.
  - b) Use mini-gel apparatus (~50 ml total gel volume).
  - c) Prepare samples by mixing together the following:
    - 2  $\mu$ l DNA sample
    - 1  $\mu$ l Orange G
    - 7  $\mu$ l ddH<sub>2</sub>O
  - d) Use 2  $\mu$ l 1 Kb Plus DNA ladder.
3. If digested completely, run samples on medium-sized gel (~100 mL total gel volume).
  - a) Use 1% agarose gel with 5  $\mu$ l Ethidium Bromide / 100  $\mu$ l agarose gel.
  - b) Add 1  $\mu$ l Orange G / 10  $\mu$ l digested sample.
  - c) Use 5  $\mu$ l 1 Kb Plus DNA ladder.
  - d) Run gel for ~4 hr at 100 V.
4. Excise 1.3 Kb band from gel (2<sup>nd</sup> band from bottom). Take pictures of before & after.
5. Extract DNA from Gel using Qiagen QIAEX II Agarose Gel Extraction Protocol.
6. Quantify sample using lambda ladder.
  - a) Use 1% agarose gel with 5  $\mu$ l Ethidium Bromide / 100  $\mu$ l agarose gel.
  - b) Use mini-gel apparatus (~50 ml total gel volume).
  - c) Prepare samples by mixing together the following:
    - 5  $\mu$ l DNA sample
    - 2  $\mu$ l Orange G
    - 3  $\mu$ l ddH<sub>2</sub>O
  - d) Prepare  $\lambda$  ladder by mixing together the following and incubating for 10 minutes at 65°C:
    - 1  $\mu$ l  $\lambda$  marker
    - 2  $\mu$ l Orange G
    - 7  $\mu$ l ddH<sub>2</sub>O
  - e) Use 5  $\mu$ l  $\lambda$  ladder.
7. Compare sample band with  $\lambda$  ladder bands to determine concentration.