

Immunoprecipitation Protocol

1. Lyse Cells

Prepare Lysis Buffer (50ml): 10% Glycerol, 50mM Tris pH 7.7, 150 mM NaCl, 0.5% NP40, 1X DTT, 1X Protease Inhibitors, 0.1 μ M Okadaic Acid or Microcysteine. Keep everything on ice.

- a. Add appropriate volume of lysis buffer to cell pellet
- b. Use 22 gauge needle and syringe to lyse cells
- c. Add sample of lysate to SDS for loading control
- d. If volume <1.5 ml, spin top speed in cold room centrifuge for 30 minutes.

2. Prepare Beads

Depending upon antibody and beads wanted, wash 5 μ l of solid beads 1X in lysis buffer. Pre-coupled antibodies work better for western blotting after IP.

3. Add supernatant from centrifugation to beads (remember all controls) and incubate beads with supernatant in cold room, rotating, for 2 hours.

4. Wash beads 4X with lysis buffer. Careful not to lose beads during washes.

5. Prepare SDS-Page gel for western blotting, or proceed to other assays.